

Retractors



WEST RETRACTOR



PN0001



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0001	West Retractor	135	18	3 x 4 Sharp

WEISLANDER RETRACTOR



PN0014



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0014	Weislander Retractor	135	18	3 x 4 Blunt



PN0951



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN1094	Weitlaner Retractor	110	15	2 x 3 Sharp
PN0951	Weitlaner Retractor	110	15	2 x 3 Semi-Sharp
PN0952	Weitlaner Retractor	110	15	2 x 3 Blunt

TRIVERS RETRACTOR



PN0015



PN0016



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0015	Travers Retractor	210	26	4 x 5 Blunt
PN0016	Travers Retractor, Curved	200	26	4 x 5 Blunt



PN0012



PN0009




CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0012	Mollison Retractor, Child	130	12	2 x 2 Sharp
PN0009	Mollison Retractor, Adult	170	15	4 x 4 Sharp

NORFOLK & NORWICH RETRACTOR




PN0007



1/1



PN0021



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0007	Norfolk & Norwich Retractor	220	50	4 x 5 Blunt
PN0021	Norfolk & Norwich Retractor, Heavy Rack	220	50	4 x 5 Blunt



PN0918



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0918	Norfolk & Norwich Retractor, Long Teeth	220	78	4 x 5 Blunt

JOLLS THYROID RETRACTOR



PN0212

CODE	DESCRIPTION	LENGTH
PN0212	Jolls Thyroid Retractor	155

HENDERSON RETRACTOR



PN0013



1/1



PN0989



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0013	Henderson Retractor	175	20	4 x 4 Sharp
PN0989	Henderson Retractor	175	20	4 x 4 Blunt



PN0008



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0008	Jefferson Retractor	130	16	3 x 4 Blunt

IRWIN RETRACTOR



PN0018



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0018	Irwin Retractor	175	20	4 x 4 Sharp

PLESTER JANSEN RETRACTOR



PN0443



PN0444



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0443	Plester Jansen Retractor, Left Solid Blade	135	Left 28, Right 20	Solid x 3 Blunt
PN0444	Plester Jansen Retractor, Right Solid Blade	135	Left 20, Right 28	3 Blunt x Solid



PN0445



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0445	Plester Jansen Retractor	135	20	3 x 3 Blunt

WALLACE MODIFIED PLESTER JANSEN RETRACTOR



PN0410



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0410	Wallace Modified Plester Jansen Retractor	135	20	2 x 2 Semi-Sharp

CAWTHORNE RETRACTOR



PN0017



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0017	Cawthorne Retractor, Long Posterior Tooth	170	15	4 x 4 Sharp

MAWSON RETRACTOR




PN0019



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0019	Mawson Retractor	140	12	3 x 3 Blunt



PN0557




1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0557	Juniper Retractor	150	30	Blunt

JACKSON BURROWS RETRACTOR



PN0004





1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0004	Jackson Burrows Retractor	180	20	Blunt

GELPI RETRACTOR



PN0005




PN0110




CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0005	Gelpi Retractor, Adult	185	20	Sharp
PN0110	Gelpi Spinal Retractor	280	78	Sharp




PN0067



PN0547 1/1



PN0067 1/1



PN0038 1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0547	Gelpi Retractor	185	30	Sharp
PN0067	Gelpi Retractor	195	50	Sharp
PN0038	Gelpi Retractor	200	70	Sharp

HARDMAN RETRACTOR



PN0023



1/1



PN0125




1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0023	Hardman Retractor, Equal Prongs	195	10	Sharp
PN0125	Hardman Retractor, Equal Prongs	215	50	Sharp



PN0022



PN0549



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0022	Hardman Retractor, Long Prong Left	215	Left 50, Right 30	Sharp
PN0549	Hardman Retractor, Long Prong Right	215	Left 30, Right 50	Sharp

CONE HEMI-LAMINECTOMY RETRACTOR



PN0447




PN0448



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0447	Cone Hemi-Laminectomy Retractor, Left	165	Left 17	Sharp
PN0448	Cone Hemi-Laminectomy Retractor, Right	165	Right 17	Sharp



PN0449



1/1



PN0450



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0449	Cone Hemi-Laminectomy Retractor, Left	255	Left 24	Sharp
PN0450	Cone Hemi-Laminectomy Retractor, Right	255	Right 24	Sharp

CONE RETRACTOR



PN0495



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0495	Cone Retractor	165	17	3 x 4 Sharp

CONE LAMINECTOMY RETRACTOR



PN0446




CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0446	Cone Laminectomy Retractor	265	24	3 x 4 Sharp

HARVEY JACKSON LAMINECTOMY RETRACTOR



PN0010



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0010	Harvey Jackson Laminectomy Retractor	305	45	3 x 3 Sharp

MAYO ADSON LAMINECTOMY RETRACTOR



PN0006




1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0006	Mayo Adson Laminectomy Retractor	305	30	4 x 4 Blunt

BROWSE RETRACTOR




PN0211



1/1



PN0440



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0211	Browse Retractor, Long Prong Left	195	Left 50, Right 25	4 x 4 Blunt
PN0440	Browse Retractor, Long Prong Right	195	Left 25, Right 50	4 x 4 Blunt



PN0441



PN0020



CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0441	Talbot Orbital Retractor, Left	165	Left 40 & 52, Right 13	Interchangeable x 4 Sharp
PN0020	Talbot Orbital Retractor, Right	165	Left 13, Right 40 & 52	4 Sharp x Interchangeable

CAIRNS SCALP RETRACTOR



PN0115



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0115	Cairns Scalp Retractor	140	15	3 x 3 Sharp

R.N.T.N.E.H. RETRACTOR



PN0011



1/1

CODE	DESCRIPTION	LENGTH	BLADE DEPTH	TEETH
PN0011	R.N.T.N.E.H. Retractor	170	15	3 x 4 Sharp



PN0324



1/1

Adult Blades



1/1

Sphincterotomy Blades



1/1

Centre Blade

CODE	DESCRIPTION	LENGTH
PN0324	Park Anal Retractor, with Adult, Sphincterotomy and Centre Blades	145



PN0558



PN0573

CODE	DESCRIPTION	LENGTH	BLADE DEPTH
PN0558	Eisenhammer Retractor	180	80
PN0573	Eisenhammer Retractor with Ratchet	180	80



PN0131

CODE	DESCRIPTION	LENGTH	BLADE DEPTH
PN0131	Farrands Retractor	180	140

SELLORS RETRACTOR



PN0041

CODE	DESCRIPTION
PN0041	Sellors Retractor



PN1058



PN1059

CODE	DESCRIPTION	LENGTH
PN1058	Harefield Sternal Blades	78 x 30
PN1059	Harefield Sternal Blades	120 x 30

IMA RETRACTOR



PN0305

CODE	DESCRIPTION
PN0305	IMA Retractor

PHLEBECTOMY HOOK



PN0801

CODE	DESCRIPTION
PN0801	Phlebectomy Hook No. 1
PN0802	Phlebectomy Hook No. 2
PN0803	Phlebectomy Hook No. 3